

# METRICS

	Ø mm	G	↔mm	PORT		n / f	REF	Σ		
				n.	mm <sup>2</sup>					
<b>INFILTRATION</b>										
Skatter	1.5	17	100	5	0.7	∩	DOIN1510-2	10		
	2.0	14	150	5	1.2	∩	DOIN2015-2	10		
	2.5	13	250	5	1.4	∩	DOIN2525-2	10		
<b>LIPOSUCTION</b>										
Dolfin	3.0	11	150	3	-14	∩	DOLF3150-2	10		
	3.0	11	250	3	-14	∩	DOLF3250-2	10		
	4.0	8	150	3	-30	∩	DOLF4150-2	10		
	4.0	8	250	3	-30	∩	DOLF4250-2	10		
	4.0	8	350	3	-30	∩	DOLF4350-2	10		
	5.0	6	250	3	-50	∩	DOLF5250-2	10		
	5.0	6	350	3	-50	∩	DOLF5350-2	10		
Dolfin Plus	4.0	8	250	7	-40	∩	DOPL4250-2	10		
Dolfin Pro	4.0	8	250	8	-50	∩	DOLP4250-2	10		
	4.0	8	350	8	-50	∩	DOLP4350-2	10		
	5.0	6	250	8	-70	∩	DOLP5250-2	10		
	5.0	6	350	8	-70	∩	DOLP5350-2	10		
<b>FAT GRAFTING</b>										
<b>HARVESTING</b>										
Picking	MEGA	3.5	10	180	5	6	∩	PICK3518-2	10	
	Picking II	MEGA	4.0	8	300	8	6	∩	PICK352A-2	10
	MACRO	3.5	10	230	18	1.7	∩	PICK352L-2	10	
<b>INJECTION</b>										
Fillin Rigid	MEGA	3.5	10	250	1	5	∩	FIL1270S-2	10	
	MEGA	2.5	12	200	1	3	∩	FILU2515-2	10	
	MEGA	2.7	12	200	1	3	∩	FIL2590S-2	10	
	MEGA	2.7	12	150	1	3	∩	F9027150-2	10	
	MACRO	2.0	14	150	1	1.9	∩	FIL25C5S-2	10	
	MACRO	2.0	14	150	1	1.9	∩	FIL1670S-2	10	
	MACRO	2.0	14	100	1	1.9	∩	F9020100-2	10	
	MACRO	1.5	17	100	1	1.0	∩	F9015100-2	10	
	MACRO	1.5	17	70	1	1.0	∩	F9015070-2	10	
	MICRO	1.2	18	100	1	0.9	∩	F9012100-2	10	
	MICRO	1.2	18	70	1	0.9	∩	F9012070-2	10	
	Fillin Flexible	MICRO	1.0	19	100	1	0.4	∩	F9010100-2	10
		MICRO	1.0	19	70	1	0.4	∩	F9010070-2	10
MICRO		0.8	21	80	1	0.2	∩	F9008080-2	10	
MICRO		0.8	21	60	1	0.2	∩	F9008060-2	10	
MICRO		0.8	21	40	1	0.2	∩	F9008040-2	10	
NANO		0.6	23	50	1	0.14	∩	F9006050-2	10	
NANO		0.5	25	40	1	0.1	∩	F9005040-2	10	
NANO	0.5	25	25	1	0.1	∩	F9005025-2	10		
<b>KITS</b>										
		A	B1	B2	B3	C1	C2			
Macro St'rim	Ø ↔mm	3.5 mm 220	12G 150	14G 150	17G 100	17G 70		PICS3523-2	5	
St'rim	Ø ↔mm	2.0 mm 130	19G 50	21G 60	21G 40	14G 40 21G 40		PICS2130-2	5	
Mini St'rim	Ø ↔mm	2.0 mm 90	19G 50	21G 60	21G 40	14G 40 21G 40		PICS2090-2	5	
Mini St'rim	Ø ↔mm	2.0 mm 70	19G 50	21G 60	21G 40	14G 40 21G 40		PICS2070-2	5	
<b>ACCESSORIES</b>										
Doto	DETACHABLE HANDLE FOR MACRO ST'RIM							DOTO0155-2	10	

∩ EFFICIENCY


LOW  
MEDIUM  
HIGH  
VERY HIGH

f EXTRUSION FORCE

LOW  
VERY LOW

# Be eSTRIM

PERFORM EXTREMELY SAFE & EFFICIENT FAT GRAFTING



TYPE OF FAT GRAFTING	SIZING HARVESTING CANNULAS Picking's			EQUALIZING INJECTION CANNULAS Fillin			KITS HARVESTING & INJECTION
	Ref.	Lobule Size	∩	Ø mm	f	Ref.	
MEGAFAT	Picking		HIGH	10 G	3,5	VERY LOW	—
	Picking II		VERY HIGH	12 G	2,7	LOW	
MACROFAT	Picking II		HIGH	12G	2,7	VERY LOW	Macro St'rim
				14G	2,0	VERY LOW	
				17G	1,5	LOW	
MICROFAT	—		—	18G	1,2	VERY LOW	St'rim Mini St'rim
				19G	1,0	VERY LOW	
				21G	0,8	LOW	
NANOFAT MAT <sup>1</sup> /tSVF <sup>2</sup>	—		—	23G	0,6	VERY LOW	—
				25G	0,5	LOW	

LIPU\_FG\_PORTFOLIO.CPCEN.C4(SEPTEMBER22)

The spheres are calibrated and not deformable and are only intended to illustrate a hypothesis of an adipocyte lobule size based on data of the technical files of STRIM HEALTH CARE.

- (1) Microfragmentated Adipose Tissue
- (2) Stromal Vascular Fraction (tissue)

- ∩ Suction efficiency
- f FILLIN'S extrusion forces
- Σ Total suction surface

Distributed by

**esthetic med**  
arzt- und spitalbedarf

estheticmed.ch | info@estheticmed.ch  
Tel. +41 (0)44 802 33 88 | Fax. +41 (0)44 802 33 89  
Chli Ebnet 3, 6403 Küsnacht am Rigi -Switzerland

FIND OUT MORE ABOUT IN  
[WWW.STRIMHEALTHCARE.COM](http://WWW.STRIMHEALTHCARE.COM)

**strim**  
HEALTH CARE


welcome@strimhc.com | +33 450 72 17 18 | Route des Jeunes 9, 1227 Les Acacias, Geneva -Switzerland

©2022

**strim**  
HEALTH CARE

## LIPOSUCTION & FAT GRAFTING


COMPACT PRODUCT CATALOG




# INFILTRATION

## SKATTER'S


Multi-spiraled cannulas for tumescent fluid infiltration prior to fat harvesting



Length: 100, 150 and 250 mm



Vortex Effect



360° Distribution


# LIPOSUCTION

## DOLFIN'S


Body reshaping cannulas with sharp edges ports "guillotine" design to cut adipose tissue




### DOLFIN



MEDIUM  $\eta$




260° SUCTION PRECISION




### DOLFIN PLUS



HIGH  $\eta$




260° SUCTION PRECISION




### DOLFIN PRO



VERY HIGH  $\eta$




360° FULL SUCTION



# FAT GRAFTING


## MEGAFAT



Ø mm	4.0	3.5	10G	12G	
← (mm)	300	180	250	200	150
PORTS	8	5	INJECTION		
$\eta$	VERY HIGH	HIGH			
SUCTION	260°	360°			

HARVESTING HARVESTING


## MACROFAT



Ø mm	3.5	14G	17G			
← (mm)	230	150	150	100	100	70
PORTS	18	INJECTION				
$\eta$	HIGH					
SUCTION	260°					


HARVESTING

## MICROFAT



← (mm)	100	70	100	70	80	60	40
--------	-----	----	-----	----	----	----	----

## NANOFAT



← (mm)	50	40	25
--------	----	----	----

ALL IN ONE

KITS

## macro st'rim®

MACROFAT GRAFTING KIT




### OPTIONAL:

- Connector FLL/FLL
- DOTO handle
- SKATTER

## st'rim®

MICROFAT GRAFTING KIT



## mini st'rim®

Similar kit to st'rim but with shorter harvesting cannulas - see metrics tab.